

AUTOCLAVABLE BIO-HAZARD BAGS

Ideal for sterilizing used petri dishes from coliform testing prior to disposal. Bags are bright orange-red high density polyethylene with the biohazard symbol. Bags measure 8 1/2" x 11" and have one long flap for ease in opening. Bags have a sterilization indicator on the outside.

BB-333 15.45/pk of 100

BAG HOLDER FOR BB-333 BIO-HAZARD BAGS

BB-333H 23.30/ea



AUTOCLAVABLE CLEAR BAGS

No printing on these bags. Measure 8" x 12" x 2 mil thick.

BB-812 41.60/pk of 100



SAMPLE SCOOP FOR WHIRL-PAK BAGS

Designed for collecting liquid samples. Two clips on the end will hold a variety of bag sizes. 18" handle is coated with a powder finish to provide durability.

WS-100 12.10/ea

BASKETS

Epoxy-coated, steel wire baskets hold up under the most severe laboratory conditions. They will resist rust, corrosion, solvents, acids, and alkalis. Handy for transporting and batch dipping glassware.



Size	Bottom Mesh	NCL #	Price/ea
11" x 6" x 3 1/2"	1/4" x 1"	RB-765	30.00
14" x 11" x 3 1/2"	1/4" x 1"	RB-766	58.50
20" x 11" x 3 1/2"	1/4" x 1"	RB-767	67.35



WHIRL-PAK BAGS: STERILE SAMPLE BAGS

Capacity (oz)	Size	NCL #	Price/pk of 500
4	3" x 7"	B-206	41.00
18	4.5" x 9"	B-218	57.80
24	6" x 9"	B-224	75.50
27	5" x 12"	B-227	78.90
42	6" x 15"	B-242	107.70

WRITE-ON STERILE SAMPLE BAGS, WHIRL-PAK

The white, write-on strip makes sample identification easier and more convenient. White strip may be written on with a ballpoint pen, pencil, or felt-tip marker.

Capacity (oz)	Size	NCL #	Price/pk of 500
4	3" x 7"	BN-206	43.40
18	4.5" x 9"	BN-218	61.15
24	6" x 9"	BN-224	80.60



THIO-BAGS

Whirl-pak bags with a Sodium Thiosulfate tablet in each bag. Sodium Thiosulfate tablet neutralizes chlorine so that sample is correct for coliform testing.

Capacity (ml)	Size	NCL #	Price/pk of 100
100	3" x 7"	B-100	19.30
300	4.5" x 9"	B-103	43.30

STERILE STAND-UP SAMPLE BAGS, WHIRL-PAK

Whirl-pak bags with a flat bottom to ensure stability and wire tabs hold the top open for easy access. Bags include a write-on strip for sample identification. The BS-100 contains a Sodium Thiosulfate tablet.

Capacity (oz)	Size	NCL #	Price/pk
4	3" x 7"	BS-100	22.90/100
4	3" x 7"	BS-206	56.00/500
18	4.5" x 9"	BS-218	82.30/500

