

CORNING GLASS BEAKERS

Corning 1000 series. Reinforced rims and uniform wall thickness offers optimum balance between thermal shock resistance and mechanical strength. The 250 through 4000 ml sizes have double graduations to indicate approximate contents. The 10 ml size is not graduated.

NOTE: Beakers are NOT suitable for high-accuracy measurements.

Size (ml)	Grad Range (ml)	Approximate Diameter x Height	NCL #	Price/ea	Price/pk	Price/cs
10		25 x 32 mm	BC-399	4.70	46.55/12	176.45/48
30	5 – 25	35 x 53 mm	BC-399A	4.00	39.90/12	151.20/48
50	10 – 40	42 x 56 mm	BC-400	4.00	39.90/12	151.20/48
100	20 – 80	50 x 72 mm	BC-401	4.25	42.20/12	160.00/48
150	20 – 140	57 x 84 mm	BC-402	3.85	38.30/12	145.10/48
250	25 – 200	68 x 90 mm	BC-403	4.10	40.45/12	153.25/48
400	25 – 325	77 x 110 mm	BC-404	4.80	38.05/12	180.10/48
600	50 – 500	90 x 124 mm	BC-405	6.10	30.15/6	171.45/36
1000	50 – 1000	108 x 158 mm	BC-407	11.45	56.85/6	215.40/24
1500	200 – 1400	121 x 166 mm	BC-408	13.90	45.95/4	174.05/16
2000	200 – 2000	131 x 193 mm	BC-409	22.75	75.15/4	142.40/8
3000	250 – 2500	146 x 216 mm	BC-409A	39.10	-----	183.65/6
4000	500 – 3500	160 x 250 mm	BC-410	72.10	-----	338.50/6



KIMAX GLASS BEAKERS

Reinforced rims and uniform wall thickness. Kimble 14000 series.

Size (ml)	Grad Range (ml)	Approximate Diameter x Height	NCL #	Price/ea	Price/pk	Price/cs
10		25 x 32 mm	BK-399	4.60	45.70/12	173.30/48
30	5 – 25	35 x 53 mm	BK-399A	4.00	39.90/12	151.15/48
50	10 – 40	41 x 53 mm	BK-400	4.00	39.90/12	151.15/48
100	20 – 80	50 x 71 mm	BK-401	4.20	41.50/12	157.25/48
150	20 – 140	57 x 79 mm	BK-402	4.00	39.75/12	150.55/48
250	25 – 200	68 x 88 mm	BK-403	4.00	39.75/12	150.55/48
400	25 – 325	77 x 106 mm	BK-404	4.85	48.05/12	182.10/48
600	50 – 500	88 x 122 mm	BK-405	6.00	29.75/6	169.15/36
1000	100 – 1000	107 x 145 mm	BK-407	11.30	27.40/6	211.80/24
1500	200 – 1400	120 x 165 mm	BK-408	13.65	45.15/4	171.05/16
2000	200 – 2000	130 x 190 mm	BK-409	22.35	73.90/4	140.00/8
4000	500 – 3500	161 x 250 mm	BK-410	70.80	-----	332.45/6



GLASS BEAKER ASSORTMENT

Consists of one each of the following sizes: 50 ml, 100 ml, 250 ml, 600 ml and 1 liter.

BC-499

28.30/pk



HEAVY-DUTY GLASS BEAKERS, WITH SPOUT.

Corning 1003 series. Heavier walls than the regular beakers, specifically designed to offer the best mechanical strength under harsh conditions.

Size (ml)	Grad Range (ml)	Approximate Diameter x Height	NCL #	Price/ea	Price/pk	Price/cs
150	20 – 140	57 x 89 mm	BC-420	7.70	76.25/12	288.95/48
250	25 – 200	68 x 90 mm	BC-421	7.60	75.30/12	285.35/48
400	25 – 325	77 x 110 mm	BC-422	9.90	98.15/12	372.05/48
600	50 – 500	90 x 124 mm	BC-423	11.55	57.20/6	325.05/36
1000	50 – 1000	108 x 156 mm	BC-424	22.05	109.25/6	413.95/24
2000	200 – 2000	131 x 180 mm	BC-425	45.00	148.75/4	281.80/8
4000	500 – 3500	160 x 250 mm	BC-426	87.90	290.45/4	-----

HEAVY-DUTY GLASS BEAKERS, WITH SPOUT.

Kimble 14005 series. Comparable to the Corning 1003 series.

Size (ml)	Grad Range (ml)	Approximate Diameter x Height	NCL #	Price/ea	Price/pk	Price/cs
250	25 – 200	67 x 92 mm	BK-421	7.45	74.00/12	280.35/48
400	25 – 325	77 x 114 mm	BK-422	9.75	99.35/12	365.70/48
600	50 – 500	87 x 122 mm	BK-423	11.35	56.25/6	319.65/36
1000	50 – 1000	107 x 152 mm	BK-424	21.95	108.70/6	412.00/24
2000	200 – 2000	130 x 190 mm	BK-425	44.25	146.20/4	277.00/8
4000	500 – 3500	161 x 252 mm	BK-426	86.40	285.50/4	-----

TALL-FORM GLASS BEAKERS

CORNING Tall-Form Beakers, with spout.

Size (ml)	Grad Range (ml)	Approximate Diameter x Height	NCL #	Price/ea	Price/pk	Price/cs
100	20 – 80	50 x 79 mm	BC-330T	5.00	49.50/12	187.55/48
200	25 – 150	56 x 102 mm	BC-350T	5.30	52.70/12	199.70/48
300	25 – 250	64 x 118 mm	BC-360T	5.90	58.40/12	221.35/48
400	25 – 325	69 x 127 mm	BC-370T	6.40	31.75/6	180.35/36
500	50 – 450	75 x 136 mm	BC-380T	7.05	34.85/6	165.15/30
600	50 – 550	79 x 152 mm	BC-390T	9.10	45.15/6	171.05/24
1000	50 – 950	89 x 187 mm	BC-410T	12.20	60.50/6	172.05/18





DISPOSABLE POLYSTYRENE BEAKERS

50 ml pour spout, no graduations. Ideal for pipetting vessel to avoid pipetting out of a reagent bottle.

BE-50D **29.00/pk of 100** **87.50/cs of 500**

GLASS BEAKERS, WITH SPOUTS & HANDLES		
Made from Berzelius glass with solid rod handles. Both the mug and the pitcher are graduated and can withstand thermal shock.		
Size (ml)	NCL #	Price/ea
500	BK-500M	23.90
2000	BK-500P	89.50



STAINLESS STEEL BEAKERS WITH SPOUTS (NO HANDLES)		
Size (ml)	NCL #	Price/ea
600	B-520	63.75
1200	B-521	96.00
2000	B-522	109.50

BEAKER COVERS		
"WATCH GLASSES" with fire-polished edges. All are heavy duty except the 8" size.		
Diameter	NCL #	Price/ea
2.5"	WC-202	3.65
3.0"	WC-203	3.85
3.5"	WC-203A	4.35
4.0"	WC-204	3.55
5.0"	WC-205	5.00
6.0"	WC-206	6.90
8.0"	WK-208	14.85



WHITE PLASTIC BEAKERS

Low-Form Beakers, polypropylene. For general laboratory use. Superior chemical resistance. All sizes have easy-to-read silk-screened graduations, plus catalog number, size, code, resin code, maximum use temperature and "no flame" symbol. With ribs for easy stacking. Not for use on a hot plate. Autoclavable.

Size (ml)	NCL #	Price/ea	Price/pk
50	BN-500	2.05	21.40/12
100	BN-501	2.50	24.80/12
150	BN-502	3.20	33.60/12
250	BN-503	3.30	18.00/6
400	BN-504	4.35	23.60/6
600	BN-505	4.95	27.00/6
1000	BN-506	6.80	37.00/6
2000	BN-507	13.70	74.60/6
4000	BN-508	52.90	172.26/4



PLASTIC BEAKERS WITH HANDLES

Metric graduations are permanently molded on the outside.



Size (ml)	NCL #	Price/ea
250	B-508	3.90
500	B-509	4.25
1000	B-510	6.00
2000	B-511	9.05
3000	B-512	12.00
5000	B-514	26.00

CLEAR PLASTIC BEAKERS

As transparent as glass. Polymethylpentene. Nalgene. Autoclavable.

Size (ml)	NCL #	Price/ea	Price/pk	Price/cs
50	BN-510	4.20	42.35/12	120.40/36
100	BN-511	4.95	51.35/12	145.90/36
150	BN-512	6.50	67.40/12	191.50/36
250	BN-513	8.50	43.90/6	166.30/24
400	BN-514	12.60	65.35/6	247.70/24
600	BN-515	15.45	53.35/4	151.65/12
1000	BN-516	18.70	48.65/3	184.40/12
2000	BN-517	36.70	-----	132.12/4
4000	BN-518	57.75	-----	207.99/4

