



### BOD BOTTLE BRUSH

With plastic handle. Bristles are held by twisted wire which is mounted on a plastic handle. Brush part is curved and handle is offset to allow you to reach all parts of a 300 ml BOD bottle. The best thing available that we know of for cleaning BOD bottles.

**B-901 5.25/ea**

### BEAKER BRUSH


With wooden handle. For cleaning beakers, jars, wide-mouth bottles, etc. Black bristles on twisted wire, mounted on a hardwood handle. Brush size: 6" long x 2<sup>3</sup>/<sub>8</sub>" diam. Overall length: 16"

**B-900 13.50/ea**



### NAIL BRUSHES

Black nylon bristles on a plastic block. Brush size: 1<sup>1</sup>/<sub>8</sub>" x 3<sup>3</sup>/<sub>4</sub>" x 5<sup>5</sup>/<sub>8</sub>". The B-370 is autoclavable.



NCL #	Price/ea
<b>B-370</b>	<b>7.75</b>
<b>B-371</b>	<b>4.40</b>




### IMHOFF CONE BRUSH

Extra heavy tip extends 1<sup>3</sup>/<sub>4</sub>". Brush portion is 9<sup>3</sup>/<sub>4</sub>" long. First 4" is horsehair, remaining 5<sup>3</sup>/<sub>4</sub>" is brown fiber on galvanized wire. Diameter of brush is 3<sup>3</sup>/<sub>4</sub>" at tip, tapering to 4" at the base. Hardwood handle is 15" x 3<sup>3</sup>/<sub>4</sub>" diameter.

**B-980 19.00/ea**

### TEST TUBE BRUSHES

Ideal for cleaning test tubes, cylinders, and similar objects. Stiff white polyester fibers are twisted in stainless steel wire. Radial tufted end. Looped handle for hanging.



Brush Size	Overall Length	NCL #	Price/ea
1/2" x 3"	8"	<b>B-913A</b>	<b>3.75</b>
3/4" x 3 1/2"	9"	<b>B-913</b>	<b>3.75</b>
1 3/8" x 4"	10 1/2"	<b>B-914</b>	<b>4.25</b>



### LAB GLASSWARE WASHING KIT

Includes Brush Organizer (Wall Mounting) and one each of B-875, B-876, and B-878.

**B-880B 25.95/ea**



### NON-ABRASIVE LAB BRUSH

Excellent for carboys, trays, circulating baths, and flat surfaces. Will not be affected by detergents and chemicals. Brush size 2.75" W x 4.75" H x 1.5" thick. Overall length: 14 1/4"

**B-253 7.95/ea**



### TEST TUBE BRUSH

Expandable nylon twisted-in-wire bristles. Fits tube sizes 3/4" ID and larger. Overall length: 8 1/2"

**B-877 4.95/ea**

### LONG-LIFE FOAM BRUSHES

Manufacturer claims that these last 10 times longer than the average brush. We don't know about that, but they are much more durable compared to other brushes.

- ✓ New special formula foam does not absorb water and does not become soggy and limp.
- ✓ Will not scratch, chip or break laboratory glassware and will not mar plasticware.
- ✓ Long-life – under normal usage will last 10 times longer than the best laboratory brush and 20 times longer than the best sponge.

#### CYLINDER BRUSH

Overall length: 18"

**B-875 7.95/ea**

#### BEAKER BRUSH


Overall length: 9"

**B-876 6.95/ea**

#### FLEXIBLE BRUSH

For narrow-mouth flasks and other glassware. Overall length: 16"

**B-878 7.95/ea**





**VOLUMETRIC FLASK BRUSHES**

Identical in construction to the BOD bottle brush mentioned previously. For cleaning the sizes indicated.

For Cleaning	NCL #	Price/ea
200/250 ml flasks	<b>B-920</b>	<b>5.25</b>
500 ml flasks	<b>B-921</b>	<b>5.50</b>
1000 ml flasks	<b>B-922</b>	<b>5.80</b>



**FLASK BRUSH, METAL HANDLE**

White bristles.  
Brush size: 5<sup>5</sup>/<sub>8</sub>" long x 2" diam.  
Overall length: 21"

**B-911 8.50/ea**



**VOLUMETRIC FLASK BRUSH**

For 50 ml and 100 ml volumetric flasks. Soft, black bristles are twisted into galvanized wire.


**B-919 4.50/ea**



**PIPET BRUSH**

Back by popular demand. Works marginally on 20 ml and 25 ml volumetric pipets. Does not work on other sizes or on serological pipets. Totally worthless. We offer them because everyone else does.

**B-925 3.50/ea**



**CYLINDER BRUSHES**

Nylon bristles with twisted wire handle.

Brush Size	Overall Length	NCL #	Price/ea
2" x 5"	17"	<b>B-915</b>	<b>4.45</b>
2 <sup>3</sup> / <sub>8</sub> " x 6"	20"	<b>B-916</b>	<b>9.75</b>



**CARBOY BRUSH**

For cleaning the insides of carboys and water bottles. Black nylon bristles in twisted wire, mounted on a metal swivel shank, attached to a 12" long wood handle. Swivel shank allows you to reach all parts of bottle interior. Overall length is 23" fully extended.


**B-912 10.65/ea**



**SMALL BOTTLE BRUSH**


A smaller version of the Carboy Brush. Brush is 1<sup>3</sup>/<sub>8</sub>" x 4" long, connected by a swivel head to a 10<sup>1</sup>/<sub>2</sub>" long handle. Overall length is 12".

**B-912A 5.00/ea**



**BABCOCK BOTTLE BRUSH**


**B-926      3.30/ea**



**LARGE BRUSH**

3" diameter. Extremely stiff nylon bristles with a heavy duty 4 ft twisted wire handle.


**B-930      19.95/ea**



**BURET BRUSHES**

With twisted wire handle.

For Cleaning	NCL #	Price/ea
25 ml burets	<b>B-905</b>	<b>3.50</b>
50 ml burets	<b>B-906</b>	<b>3.50</b>



**BRUSH HOLDER**


High-impact 18" plastic channel with 6 hooks. Extends brush life, cleaning performance, and best of all, keeps brushes easy to find. Rust-free and easy to keep clean. Has mounting holes on 16" centers. Every lab should have at least one.

**B-989      24.00/ea**



**BALANCE DUSTING BRUSH**

**B-265**      **3.50/ea**



**CAMEL HAIR BALANCE BRUSH**

Overall length: 7". Brush length: 1".


Brush Width	NCL #	Price/ea
1/2"	<b>B-266</b>	<b>7.50</b>
1"	<b>B-267</b>	<b>10.20</b>



**ANTI-STATIC BALANCE BRUSH**

New design blends natural hair and synthetic fibers that remove static as it cleans. Because it won't scratch, it is ideal for cleaning balances, microscopes, cuvettes, electronics, and computer monitors. 1" brush with 5" composite handle.

**B-321**      **17.95/ea**



**ANTI-STATIC CLEANING BRUSHES**

Goat hair bristles. Overall length: 2<sup>5</sup>/<sub>8</sub>".

Brush Width	NCL #	Price/ea
1 <sup>1</sup> / <sub>4</sub> "	<b>B-268</b>	<b>16.30</b>
2 <sup>1</sup> / <sub>2</sub> "	<b>B-269</b>	<b>23.25</b>

**SCRUB BRUSHES**



**TABLE DUSTING BRUSH**

For dusting countertops, etc. Bristles are 2<sup>1</sup>/<sub>2</sub>" long. Overall length 12".

**B-918**      **9.75/ea**



Bristles	Brush Size	Overall Length	NCL #	Price/ea
nylon	5" x 5"	20"	<b>B-350</b>	<b>17.90</b>
tampico	4 <sup>3</sup> / <sub>4</sub> " x 4 <sup>3</sup> / <sub>4</sub> "	9"	<b>B-351</b>	<b>7.90</b>