

CLASS A GLASS BURETS WITH BUILT-IN FUNNEL TOP

With PTFE stopcock.



Size	Sub-divisions	NCL #	Price/ea
5 ml	0.01 ml	BK-950	281.75
10 ml	0.02 ml	BK-951	281.75

CLASS A GLASS BURETS WITH PTFE STOPCOCK

Size	Sub-divisions	NCL #	Price/ea
10 ml	0.05 ml	BK-954	181.35
25 ml	0.1 ml	BK-955	190.05
50 ml	0.1 ml	BK-956	200.55



REPLACEMENT PTFE STOPCOCK ASSEMBLY FOR ABOVE BURETS, 2 mm BORE

BK-802 56.80/ea

SELF-ZEROING BURETS, CLASS A

Serialized and certified, with 3-way PTFE stopcock and reservoir bottle.

Size	Sub-divisions	Tolerance	Reservoir	NCL #	Price/ea
10 ml	0.05 ml	±0.02 ml	1000 ml	BK-954A	580.75
25 ml	0.10 ml	±0.03 ml	1000 ml	BK-955A	580.75
50 ml	0.10 ml	±0.05 ml	2000 ml	BK-956A	603.50



ECONOMY SELF-ZEROING BURET

Class A glass buret with polyethylene reservoir. White 0.1 ml graduations. Squeezing reservoir fills buret, and overflow drains back into reservoir. Large, stable plastic base can be detached from reservoir.

Reservoir Capacity	Buret Capacity	NCL #	Price/ea
500 ml	10 ml	BB-710	175.75
500 ml	25 ml	BB-725	204.25
1000 ml	10 ml	BB-810	181.45
1000 ml	25 ml	BB-825	209.95

PLASTIC BURET WITH PTFE STOPCOCK

Levels of colorless liquids are easily read through this crystal clear plastic buret. PTFE stopcock does not need lubrication. Not designed for exposure to alcohol, organic solvents, hydrofluoric acid, or ammonium hydroxide. Can tolerate other acids up to 5% by weight or bases up to 30% by weight. Transparent.

Size	Sub-divisions	Tolerance	NCL #	Price/ea
25 ml	0.1 ml	± 0.06 ml	BN-150	84.95
50 ml	0.1 ml	± 0.1 ml	BN-151	96.35
100 ml	0.2 ml	± 0.2 ml	BN-152	108.05



BURET ACCESSORIES



BULB FOR SELF-ZEROING BURETS

Aspirator bulb with 6" rubber tubing.
B-963 **28.20/ea**



GLASS BURET FUNNEL

FK-500 **8.10/ea**



PLASTIC BURET FUNNEL

BN-971 **2.95/ea**

BURET SUPPORT		
With weighted base, rod, and double buret clamp.		
Item	NCL #	Price/ea
complete unit	B-975	91.00
clamp only	B-976	29.00
rod & base only	B-977	68.00



SMALL ORIFICE CLEANER



Great for pipets and buret tips! A pin vise holds stainless steel wire securely while cleaning. A 7/8" diameter ring magnet on the handle keeps the cleaner handy and is useful for picking up wires. Includes wire assortment which can be stored inside the hollow handle. Length without wire is 4".

SB-800 **69.00/ea**

REPLACEMENT WIRES FOR ORIFICE CLEANER

Package of 30 stainless steel wires. 2 inches long.

SB-801 **32.40/set**



BURET CONVERSION KIT

Includes 1000 ml low-density polyethylene bottle, four feet of clear NALGENE PVC plastic tubing, polypropylene fill-tube and adapter. Makes any plastic or glass 25 or 50 ml buret self-zeroing, economically and easily. Saves you money over costly glass self-zeroing burets. No more funneling. Hazardous, wasteful spills are eliminated. Can be used with almost all aqueous titrants. Set adjustable fill-tube to zero mark on buret. Squeeze bottle to fill buret just above zero mark. Release, and excess liquid returns to bottle and titrant falls to the zero mark. To switch to a different size buret, simply reverse the adapter.

BN-970 **22.75/ea**

