



ECONOMICALLY PRICED CENTRIFUGE WITH TRUNION HEAD

Well worth your consideration if you don't have 2 grand to blow on a centrifuge. This is the only centrifuge we know of with a trunion head under \$2000 list price. This is a well-built, heavy-duty unit of stainless steel construction, without the frills of a medical centrifuge. Many operators have found this unit useful for sludge solids testing. Supplied with 6 plastic tapered centrifuge tubes which are calibrated in percent. Cannot be used with other centrifuge tubes. Tube size supplied is 19 mm.

CR-100 **875.00/ea**

Features:

- Stainless Steel Construction
- 6-Place Aluminum Trunion Head
- 15 Minute Timer
- Single Speed Drive (3000 rpm)
- Emergency Toggle Switch

EXTRA CENTRIFUGE TUBES FOR RAVEN CENTRIFUGE

Pack of 12 extra plastic tapered centrifuge tubes for CR-100 centrifuge. Calibrated in percent.

Diameter	NCL #	Price/pk
19 mm	CR-102	81.00
17 mm	CR-102A	81.00



CENTRIFUGE PACKAGE

Includes the above centrifuge, pack of 6 centrifuge tubes, three 1000 ml Raven settleometers, and one 16 ft. Raven coretaker.

CR-101 **995.00/ea**

