

THERMO SCIENTIFIC ORION TOTAL CHLORINE COMBINATION PROBE

Measures Total Chlorine, as in Wastewater. Cannot be used for Free Chlorine. No separate reference probe required. Thermo Scientific Orion states that this probe will detect concentrations of residual chlorine down to 0.01ppm. BNC connector.

O-9770BN 540.00/ea

POLISHING STRIPS

For polishing electrodes of the Total Chlorine probe.

O-948201 74.70/pk of 24

**REAGENTS FOR USE WITH CHLORINE PROBE**

Reagent	Manufactured by	NCL #	Price
Acid Reagent for Cl ₂	Thermo Scientific Orion	O-977011	69.30/475 ml HM
Acid Reagent for Cl ₂	NCL	A-13	22.90/500 ml HM
Iodide Reagent for Cl ₂	Thermo Scientific Orion	O-977010	80.10/5x60 ml
Iodide Reagent for Cl ₂	NCL	I-55	22.90/pint