



### MOHR PINCH CLAMPS

Stainless steel Mohr type pinch clamps for rubber tubing.

Maximum Tubing Size	NCL #	Price/ea
3/8" Outside Diameter	C-250	2.95
1/2" Outside Diameter	C-256	4.10
5/8" Outside Diameter	C-257	4.05



### DAY PINCH CLAMPS

White-plated, for rubber tubing up to 3/8" OD.

**C-251 2.40/ea**



### HOSE CLAMPS

Screw type hose clamp.

Maximum Tubing Size	NCL #	Price/ea
1/2" Outside Diameter	C-252	3.45
1" Outside Diameter	C-253	3.95



### HOSE CLAMP

Useful for connecting pressure or vacuum lines of rubber tubing to glass or metal tubing. Band-ribbed for greater strength. Can be tightened with a screwdriver. For tubing up to 5/8" OD.

**CL-199 1.90/ea**



### EXTENSION CLAMP

With plasticized rounded jaws, which holds objects up to 2 1/4" in diameter. Overall length 8". Requires clamp holder.

**C-261 19.80/ea**

### CLAMP HOLDER

For connecting extension clamp C-261 and similar clamps to a ring stand or support stand.

**C-260 11.50/ea**



### DOUBLE BURET CLAMP

Clamp only.

**B-976 29.00/ea**



### SUPPORT RINGS WITH CLAMP ATTACHED



Diameter	NCL #	Price/ea
2"	C-272	12.70
3"	C-273	14.10
4"	C-274	15.50
5"	C-275	17.05

## TUBING CLAMPS

Clamps will not distort, rust, or corrode. Select clamp below for desired tubing size.

Twelve-position ratchet controls flow from full open to complete closure. For tubing up to  $\frac{3}{4}$ " OD. 6 per pack.

**C-352**      **16.45/pk**



Fifteen-position adjustment controls flow from full open to complete closure. Fits tubing from  $\frac{1}{8}$ " to  $\frac{7}{16}$ " outside diameter. 12 per pack.

**C-351**      **9.65/pk**

Single position on/off clamp  
For tubing under  $\frac{3}{16}$ " outside diameter. 10 clamps per pack.

**C-350**      **2.70/pk**



## THERMOMETER HOLDER

This convenient holder clips onto the edge of a beaker or other vessel and holds a thermometer vertically near the side of the container for safe, easy reading. The thermometer is held in place by small "fingers" that will expand to hold thermometers ranging from  $\frac{1}{8}$ " to  $\frac{7}{16}$ " in diameter. Made from polypropylene, it withstands temperatures up to 121°C (250°F). Autoclavable. 3 per pack.

**TB-906**      **12.65/pk**

## TUBING HOLDER

Hold flexible tubing steady in beakers or pitchers with this polycarbonate clip. Fits securely onto the edge and has a 2" long arm with a hole to accept tubing up to  $\frac{3}{8}$ " outside diameter. 6 per pack.

**TB-808**      **15.85/pk**

