



Component	NCL #	Price
Complete Set	<b>CF-100</b>	<b>525.00/ea</b>
Pen Only	<b>CF-101</b>	<b>249.85/ea</b>
Light Box Only	<b>CF-102</b>	<b>134.35/ea</b>
Grid Only (pk of 5)	<b>CF-103</b>	<b>19.95/pk</b>
AC Adapter	<b>CF-104</b>	<b>38.85/ea</b>
Magnifier Only	<b>CF-105</b>	<b>136.70/ea</b>

### COLONY COUNTER SYSTEM

The complete set combines a pen-style colony counter, a mini-magnifier, a mini light box with an AC line current adapter, and a counting grid.

The mini light box offers color corrected fluorescent viewing surface that is ideal for viewing slides, plates, mini-gels, petri dishes or x-rays. Operates with 4 C batteries (not included) or with 6V adapter (included). A removable counting grid is included. Dimensions: 6<sup>1</sup>/<sub>4</sub>" x 5<sup>3</sup>/<sub>8</sub>" x 1<sup>3</sup>/<sub>4</sub>"

In a single motion this pen marks, counts and confirms with a beep. The lightweight 1.2 oz pen-style colony counter features a four-digit LED display that shows cumulative count. A selector switch allows for counting up or down to correct errors or to check the count by counting backwards. Pen is powered by a 3V battery. The long life felt-tip marker works on plastic or glass and is oil and water resistant but removable with ethyl alcohol. Comes with red and black markers.

Dimensions: 6" x 3<sup>3</sup>/<sub>4</sub>" x 3<sup>3</sup>/<sub>4</sub>"

This free standing magnifier allows enhanced viewing of slides, petri dishes, multi-well plates, gels or text. The precision molded lens offers 1.75x magnification.



Component	NCL #	Price
Counting Pen	<b>CC-133</b>	<b>16.95/ea</b>
Repl. Battery	<b>CC-104</b>	<b>8.25/ea</b>

### ECONOMY COUNTING PEN

Counter-Pen is a combination marker and digital counter designed for the lab and plant. Each time a black mark is made with the felt-tipped pen, it automatically sounds a beep, and LCD displays the counts. Use it to count cell growth on petri dishes, inventory parts, and tally samples received.

Touch, sound, and count occur simultaneously. Durable ABS plastic unit provides audio/visual counts to 99999 with reset button, on/off switch, and pen cap. Size is 6 x 2<sup>2</sup>/<sub>3</sub> inches and weight is 1 ounce.