

CLEAR PLASTIC GRADUATED CYLINDERS

Nalgene PMP graduated cylinders, completely transparent, with blue polypropylene base. Can be chemically sterilized without affecting accuracy. Shatterproof and as clear as glass. Easy to read. Can be autoclaved, but accuracy will be affected.



Size (ml)	Subdivisions	NCL #	Price/ea	Price/cs
10	0.2	CN-590	14.40	311.40/24
25	0.5	CN-591	18.05	265.95/18
50	1.0	CN-592	19.80	291.95/18
100	1.0	CN-593	22.50	217.70/12
250	2.0	CN-594	29.90	293.15/12
500	5.0	CN-595	39.05	255.70/8
1000	10.0	CN-596	44.80	242.30/6
2000	20.0	CN-597	74.95	269.90/4

ECONOMY GRADE CLEAR PLASTIC GRADUATED CYLINDERS

Nalgene PMP graduated cylinders, completely transparent. One-piece construction. For those of you who want CHEAP graduated cylinders, HERE THEY ARE!

Size (ml)	Subdivisions	NCL #	Price/ea	Price/cs
10	0.2	CB-590	5.15	100.15/24
25	0.5	CB-591	6.20	90.35/18
50	1.0	CB-592	7.80	112.85/18
100	1.0	CB-593	9.35	90.65/12
250	2.0	CB-594	12.45	120.15/12
500	5.0	CB-595	17.40	112.25/8
1000	10.0	CB-596	21.70	105.15/6



PLASTIC GRADUATED CYLINDERS WITH HANDLES

White polypropylene graduated cylinders have built-in handles to make pouring easier and safer. The handle provides better control when pouring and stopping. The following two sizes are available.

Size (ml)	NCL #	Price/ea
1000	CP-576	23.00
2000	CP-577	45.00

GLASS GRADUATED CYLINDERS, PYREX

Graduated glass cylinders with single white metric scale.

Corning 3024 series. Calibrated "to deliver." Class B.

Size (ml)	Grad Div	NCL #	Price/ea	Price/pk	Price/cs
10	0.2	CC-570	16.65	54.95/4	312.15/24
25	0.5	CC-571	23.45	77.55/4	330.65/18
50	1.0	CC-572	24.40	80.60/4	343.65/18
100	1.0	CC-573	29.10	96.15/4	273.30/12
250	2.0	CC-574	40.25	132.95/4	377.85/12
500	5.0	CC-575	56.25	—	352.30/8
1000	10.0	CC-576	74.00	—	259.00/4



GLASS GRADUATED CYLINDERS, KIMAX

Graduated glass cylinders with single white metric scale.

Kimble 20024 series. Calibrated "to deliver." Class B.

Size (ml)	Grad Div	NCL #	Price/ea	Price/pk	Price/cs
10	0.2	CK-570	17.10	56.45/4	160.35/12
25	0.5	CK-571	24.00	79.20/4	225.15/12
50	1.0	CK-572	24.90	82.25/4	233.80/12
100	1.0	CK-573	29.70	98.20/4	279.10/12
250	2.0	CK-574	41.10	135.80/4	386.00/12
500	5.0	CK-575	56.90	—	359.75/8
1000	10.0	CK-576	75.00	—	270.85/4



GLASS GRADUATED CYLINDER ASSORTMENT

Consists of one each of the following sizes: 10 ml, 25 ml, 50 ml, 100 ml and 250 ml.


CC-399 156.00/pk

GLASS GRADUATED CYLINDERS, LOW-FORM, RED-LINE, PYREX

Cylinders are low form for greater stability and feature a double pour spout.

Only the following two sizes are available. Class B.

Size (ml)	Grad Div	OD x Ht (mm)	NCL #	Price/ea	Price/cs
50	1.0	29 x 148	CC-572LF	39.80	373.85/12
100	2.0	39 x 157	CC-573LF	45.95	431.65/12



GLASS GRADUATED CYLINDERS, RED-LINE, KIMAX

Same as the Kimble 20024 series, except with a red stripe, which some technicians find makes it easier to read the level of liquid. Calibrated "to deliver." Class B.

Size (ml)	NCL #	Price/ea	Price/cs
10	CK-570R	22.25	209.15/12
25	CK-571R	34.40	323.25/12
50	CK-572R	32.50	352.50/12
100	CK-573R	38.25	359.35/12
250	CK-574R	58.95	553.80/12
500	CK-575R	70.15	439.25/8
1000	CK-576R	104.25	326.40/4

CLASS A GLASS GRADUATED CYLINDERS, KIMAX

Serialized and certified. Each cylinder has an individual serial number and is supplied with a Certificate of Graduation Accuracy. Calibrated "to deliver."



Size (ml)	Grad Div	Tolerance (ml)	NCL #	Price/ea	Price/cs
10	0.1	0.10	CK-570A	84.00	454.50/6
25	0.2	0.17	CK-571A	90.00	487.80/6
50	1.0	0.25	CK-572A	90.00	491.05/6
100	1.0	0.40	CK-573A	99.00	544.75/6
250	2.0	0.80	CK-574A	105.50	394.00/4
500	5.0	1.30	CK-575A	121.60	461.10/4
1000	10.0	2.50	CK-576A	149.80	283.90/2

GLASS MIXING CYLINDERS WITH GROUND GLASS STOPPER

Graduated glass cylinders with ground glass stoppers.



Size (ml)	NCL #	Price/ea	Price/cs
25	CK-651	34.25	113.25/4
50	CK-652	43.25	141.70/4
100	CK-653	49.55	163.75/4
250	CK-654	72.90	240.95/4
500	CK-655	94.35	443.15/6
1000	CK-656	108.20	410.05/4

GLASS UNGRADUATED CYLINDERS

For use with hydrometers.



Size (ml)	Height (mm)	NCL #	Price/ea
100	254	CC-680	19.60
250	340	CC-681	26.05
500	390	CC-682	39.75
1000	460	CC-683	49.65

PLASTIC GRADUATED CYLINDERS

Polypropylene, Class B, translucent, ring-molded graduations overprinted in blue. Hexagonal base.



Size (ml)	NCL #	Price/ea	Price/pack
10	CB-646	4.40	43.20/12
25	CB-647	5.40	52.75/12
50	CB-648	6.50	63.90/12
100	CB-649	8.00	78.50/12
250	CB-650	11.30	55.45/6
500	CB-651	16.45	80.65/6
1000	CB-652	20.60	110.05/6
2000	CB-653	37.00	99.90/3