



EVAPORATING DISHES Coors porcelain							
Capacity (ml)	Top (mm)	Height (mm)	Coors #	NCL #	Price/ea	Price/pk	Price/cs
35	60	24	60196	EC-498	9.05	93.80/12	355.40/48
70	75	30	60197	EC-499	9.15	94.60/12	358.55/48
80	80	30	60198	EC-500	9.85	51.15/6	193.70/24
120	94	42	60200	EC-501	12.90	66.75/6	253.00/24
150	100	42	60201	EC-502	14.85	76.90/6	218.50/18
250	115	45	60202	EC-503	18.25	94.70/6	269.10/18



SHALLOW FORM EVAPORATING DISHES Coors porcelain							
Capacity (ml)	Top (mm)	Height (mm)	Coors #	NCL #	Price/ea	Price/pk	Price/cs
30	70	16	60230	EC-530	10.60	54.85/6	207.70/24
50	80	20	60231	EC-531	11.35	58.70/6	222.55/24
75	90	22	60232	EC-532	13.90	71.95/6	204.50/18
100	105	23	60233	EC-533	14.90	77.10/6	219.20/18



CASSEROLE DISHES Coors porcelain				
Capacity (ml)	Coors #	NCL #	Price/ea	Price/cs
140	60059	CC-187	25.10	246.70/12
210	60060	CC-188	29.90	293.50/12
500	60062	CC-190	48.50	238.00/6



ALUMINUM WEIGHING DISHES

For TSS filters. Weigh filter in the aluminum weighing dish before using. Transfer dish and filter to drying oven, desiccator, and re-weigh filter in the dish. **NOT SUITABLE for use in muffle furnaces.**
100 dishes per pack, 1000 dishes per case.

Diameter	NCL #	Price/pk	Price/cs
57 mm	W-355	10.65	79.80
70 mm	W-357	14.80	122.50
125 mm	W-399	39.45	299.00