



GLASS ERLENMEYER FLASK ASSORTMENT

Consists of one each of the following sizes: 50 ml, 125 ml, 250 ml, 500 ml and 1 liter.

FC-699 40.80/pk

NARROW-MOUTH ERLENMEYER FLASKS, PYREX

Size (ml)	NCL #	Price/ea	Price/pk	Price/cs
25	FC-849	5.50	54.75/12	207.40/48
50	FC-850	5.40	53.50/12	202.75/48
125	FC-851	5.20	51.40/12	194.80/48
250	FC-852	4.95	49.15/12	186.25/48
500	FC-853	6.20	30.70/6	174.35/36
1000	FC-854	10.05	49.85/6	188.95/24
2000	FC-855	23.70	78.25/4	148.30/8




NARROW-MOUTH ERLENMEYER FLASKS, KIMAX

Size (ml)	NCL #	Price/ea	Price/pk	Price/cs
25	FK-849	4.90	48.85/12	185.05/48
50	FK-850	4.90	48.85/12	185.05/48
125	FK-851	4.70	46.65/12	176.70/48
250	FK-852	4.90	48.35/12	183.30/48
500	FK-853	6.15	30.45/6	173.15/36
1000	FK-854	10.00	49.50/6	187.50/24
2000	FK-855	23.30	76.90/4	145.70/8

WIDE-MOUTH ERLENMEYER FLASKS, PYREX

Size (ml)	NCL #	Price/ea	Price/pk	Price/cs
125	FC-861	6.50	64.70/12	245.10/48
250	FC-862	5.90	58.45/12	221.50/48
500	FC-863	7.40	36.80/6	209.15/36
1000	FC-864	12.50	62.10/6	235.25/24
2000	FC-865	25.10	82.90/4	157.05/8



WIDE-MOUTH ERLENMEYER FLASKS, KIMAX

Size (ml)	NCL #	Price/ea	Price/pk	Price/cs
125	FK-861	5.50	54.55/12	206.65/48
250	FK-862	5.85	58.05/12	219.90/48
500	FK-863	7.40	36.55/6	207.85/36
1000	FK-864	12.30	61.00/6	231.10/24
2000	FK-865	23.05	76.20/4	144.40/8



KJELDAHL FLASKS

With long neck for Kjeldahl digestions.

Size (ml)	NCL #	Price/ea	Price/pk
100	FK-111	18.25	90.45/6
500	FK-555	25.20	124.95/6
800	FK-800	22.95	113.80/6

KJELDAHL FLASKS, KOOL-GRIP®

Borosilicate flask with thicker heat insulating Kool-Grip® jacket to accommodate for the higher temperatures used in Kjeldahl digestions.

Size (ml)	NCL #	Price/ea
800	FK-800G	45.00



SUCTION FLASKS

GLASS SUCTION FLASKS, PYREX

For vacuum filtrations. Flasks take 3/8" inside diameter vacuum tubing.

Size (ml)	Stopper Size	NCL #	Price/ea	Price/pk	Price/cs
500	7 or 7 1/2	FC-800	26.55	131.70/6	374.30/18
1000	7 1/2 or 8	FC-801	53.55	265.40/6	502.85/12
2000	8 1/2	FC-802	92.00	-----	998.75/12

GLASS SUCTION FLASKS, KIMAX

For vacuum filtrations. Flasks take 3/8" inside diameter vacuum tubing.

Size (ml)	Stopper Size	NCL #	Price/ea	Price/pk
500	7 or 7 1/2	FK-800A	33.80	167.60/6
1000	7 1/2 or 8	FK-801	53.35	264.55/6
2000	8 1/2	FK-802	87.50	-----



PLASTIC SUCTION FLASKS

For vacuum filtrations. Molded in one piece to withstand full vacuum to 30 in Hg. Tubulation at 45° angle prevents tipping.

Size (ml)	Fits tubing inside diameter	NCL #	Price/ea
500	1/4" to 3/8"	FN-820	18.60
1000	1/4" to 3/8"	FN-821	29.50
1700	3/8"	FN-822	35.40



PYREX CLASS A VOLUMETRIC FLASKS WITH GROUND GLASS NECKS & HOLLOW GLASS STOPPERS

Size (ml)	Stopper Size	NCL #	Price/ea	Price/pk	Price/cs
50	9	FC-880	35.35	175.35/6	332.25/12
100	13	FC-881	37.75	187.10/6	354.50/12
200	13	FC-881A	43.60	216.15/6	409.55/12
250	16	FC-882	47.00	233.00/6	441.45/12
500	19	FC-883	58.40	289.50/6	548.50/12
1000	22	FC-884	70.75	332.20/6	-----
2000	27	FC-885	89.40	295.40/4	-----



HEAVY DUTY WIDE-MOUTH PYREX CLASS A VOLUMETRIC FLASKS WITH GROUND GLASS NECKS & HOLLOW GLASS STOPPERS



Size (ml)	Stopper Size	NCL #	Price/ea	Price/cs
10	13	FC-601	32.70	153.45/6
25	13	FC-602	32.65	153.45/6
50	13	FC-603	34.05	159.90/6
100	16	FC-604	51.00	239.50/6
200	19	FC-605	57.70	271.05/6
250	19	FC-606	64.40	302.45/6

KIMAX CLASS A VOLUMETRIC FLASKS WITH GROUND GLASS NECKS & SOLID PENNYHEAD GLASS STOPPERS

The FK-880W is a wide neck class A flask, which makes filling and pouring much easier.

Size (ml)	Stopper Size	NCL #	Price/ea	Price/pk	Price/cs
50	9	FK-880	35.80	177.45/6	336.20/12
50	13	FK-880W	35.85	168.25/6	-----
100	13	FK-881	38.20	189.35/6	358.80/12
200	16	FK-881A	44.10	218.60/6	414.20/12
250	16	FK-882	47.60	236.00/6	447.15/12
500	19	FK-883	59.10	292.80/6	554.75/12
1000	22	FK-884	71.55	336.05/6	-----
2000	27	FK-885	90.40	283.05/4	-----



PYREX CLASS A VOLUMETRIC FLASKS WITH GROUND GLASS NECKS & PLASTIC STOPPERS



Size (ml)	Stopper Size	NCL #	Price/ea	Price/pk	Price/cs
50	9	FC-890	30.85	152.90/6	289.75/12
100	13	FC-891	32.00	158.60/6	300.45/12
200	13	FC-891A	37.35	185.15/6	350.75/12
250	16	FC-892	39.65	196.65/6	372.65/12
500	19	FC-893	51.00	252.75/6	478.85/12
1000	22	FC-894	62.25	292.20/6	-----
2000	27	FC-895	78.65	259.90/4	-----

KIMAX CLASS A VOLUMETRIC FLASKS WITH GROUND GLASS NECKS & PLASTIC STOPPERS

Size (ml)	Stopper Size	NCL #	Price/ea	Price/pk	Price/cs
50	9	FK-890	31.20	154.60/6	292.90/12
100	13	FK-891	32.35	160.45/6	304.05/12
200	16	FK-891A	37.80	187.30/6	354.95/12
250	16	FK-892	40.15	198.95/6	376.95/12
500	19	FK-893	51.60	255.65/6	484.35/12
1000	22	FK-894	62.95	312.00/6	-----
2000	27	FK-895	79.55	262.85/4	-----

STOPPERS & SNAP CAPS

NOTE: Stoppers and snap caps are NOT interchangeable. Stoppers are for flasks with ground-glass joints. Snap caps will not work on these flasks. Conversely, stoppers will not work on flasks which are made for snap caps (these flasks DO NOT have ground-glass joints).

PLASTIC STOPPERS FOR VOLUMETRIC FLASKS

WATCH STOPPER SIZES! These stoppers are compatible with both Pyrex and Kimax volumetric flasks with ground-glass necks, but there is a difference in the 100 ml, 200 ml, and 250 ml Pyrex and Kimax flasks.



Stopper Size	NCL #	Price/ea
9	FC-911	3.20
13	FC-912	3.20
16	FC-913	4.20
19	FC-914	4.20
22	FC-915	5.30
27	FC-916	5.30

GLASS STOPPERS FOR VOLUMETRIC FLASKS, SOLID PENNYHEAD (KIMBLE)

WATCH STOPPER SIZES! These stoppers are compatible with both Pyrex and Kimax volumetric flasks with ground-glass necks, but there is a difference in the 100 ml, 200 ml, and 250 ml Pyrex and Kimax flasks.



Stopper Size	NCL #	Price/ea
9	FK-911	8.55
13	FK-912	8.70
16	FK-913	11.85
19	FK-914	12.75
22	FK-915	14.85
27	FK-916	17.10

GLASS STOPPERS FOR VOLUMETRIC FLASKS, HOLLOW-GROUND (PYREX)

WATCH STOPPER SIZES! These stoppers are compatible with both Pyrex and Kimax volumetric flasks with ground-glass necks, but there is a difference in the 100 ml, 200 ml, and 250 ml Pyrex and Kimax flasks.



Stopper Size	NCL #	Price/ea
9	FK-911A	9.30
13	FK-912A	11.00
16	FK-913A	11.00
19	FK-914A	14.90
22	FK-915A	15.90
27	FK-916A	19.85



VOLUMETRIC FLASKS WITH SNAP CAPS



PYREX CLASS A VOLUMETRIC FLASKS WITH SNAP-ON CAPS

Size (ml)	NCL #	Price/ea	Price/cs
25	FC-899	23.30	218.85/12
50	FC-900	23.95	224.85/12
100	FC-901	26.10	245.10/12
200	FC-901A	31.30	294.10/12
250	FC-902	33.70	316.45/12
500	FC-903	40.95	384.65/12
1000	FC-904	49.60	232.95/6
2000	FC-905	78.60	246.00/4

SNAP CAPS FOR PYREX VOLUMETRIC FLASKS

Size (ml)	NCL #	Price/ea	Price/cs
25	FC-899B	2.60	12.75/6
50	FC-900B	2.60	12.75/6
100	FC-901B	2.60	12.75/6
200	FC-906B	3.45	16.15/6
250	FC-902B	3.40	16.65/6
500	FC-903B	4.25	20.85/6
1000	FC-904B	4.30	21.10/6
2000	FC-905B	5.70	27.90/6



PLASTIC VOLUMETRIC FLASK

Nalgene polypropylene volumetric flask with screw cap. The 200 ml size has an additional mark at 203 ml so it may be utilized for Winkler testing for BOD. Autoclavable/Leakproof.



Size (ml)	NCL #	Price/ea	Price/cs
50	FN-201	21.55	232.55/12
100	FN-202	23.15	249.85/12
200	FN-203	26.00	280.80/12
250	FN-204	28.55	307.90/12
500	FN-205	34.50	248.40/8
1000	FN-206	44.10	238.25/6