

### BURET FUNNEL, GLASS

Top diameter: 35 mm; Stem OD: 6 mm; Stem length: 50 mm. Also fits volumetric flasks 50 ml and larger.

**FK-500 8.10/ea**



### GLASS FUNNELS WITH SHORT, NARROW STEM

Diam., Top (mm)	Outside Diam., Stem (mm)	Length, Stem (mm)	NCL #	Price/ea
35	6	50	<b>FK-500</b>	<b>8.10</b>
75	8	75	<b>FK-501B</b>	<b>12.60</b>
100	9	97	<b>FK-503B</b>	<b>16.30</b>

### GLASS FUNNELS WITH LONG, NARROW STEM

Diam., Top (mm)	Outside Diam., Stem (mm)	Length, Stem (mm)	NCL #	Price/ea
55	8	150	<b>FK-500B</b>	<b>8.75</b>
75	8	150	<b>FK-501</b>	<b>9.75</b>

### GLASS FUNNELS WITH SHORT, WIDE STEM ("POWDER FUNNELS")

Diam., Top (mm)	Outside Diam., Stem (mm)	Length, Stem (mm)	NCL #	Price/ea
60	13	35	<b>FK-502</b>	<b>9.00</b>
100	18	35	<b>FK-503</b>	<b>14.70</b>



### BUCHNER FUNNELS, GLAZED PORCELAIN

Diameter	NCL #	Price/ea	Price/cs
5.5 cm	<b>FC-905C</b>	<b>32.50</b>	<b>212.70/8</b>
7.0 cm	<b>FC-907</b>	<b>42.75</b>	<b>209.85/6</b>
9.0 cm	<b>FC-909</b>	<b>55.55</b>	<b>181.75/4</b>
11.0 cm	<b>FC-911A</b>	<b>69.20</b>	<b>226.45/4</b>



### BUCHNER FUNNELS, POLYPROPYLENE

Diameter	NCL #	Price/ea
5.5 cm	<b>FN-905</b>	<b>12.60</b>
7.0 cm	<b>FN-907</b>	<b>16.05</b>
9.0 cm	<b>FN-909</b>	<b>21.05</b>
11.0 cm	<b>FN-911</b>	<b>57.60</b>

## SEPARATORY FUNNELS



### PEAR-SHAPED GLASS WITH PTFE STOPCOCK & PLASTIC STOPPER

Size (ml)	NCL #	Price/ea	Price/cs
125	<b>SK-501</b>	<b>90.95</b>	<b>327.50/4</b>
250	<b>SK-502</b>	<b>108.75</b>	<b>391.50/4</b>
500	<b>SK-503</b>	<b>124.45</b>	<b>448.05/4</b>
1000	<b>SK-504</b>	<b>173.00</b>	<b>311.40/2</b>
2000	<b>SK-505</b>	<b>249.90</b>	<b>449.85/2</b>



### BURET FUNNEL, PLASTIC

Top dia: 35 mm; Stem OD: 7 mm; Stem length: 53 mm.  
Also fits volumetric flasks 50 ml and larger.

**BN-971**                      **2.95/ea**

### PLASTIC FUNNELS WITH SHORT, WIDE STEM ("POWDER FUNNELS")

For transferring powders. Parallel stem minimizes bridging of powder; external ribbing prevents air lock. Fast and efficient. Autoclavable.

Diam., Top (mm)	Outside Diam., Stem (mm)	Length, Stem (mm)	NCL #	Price/ea
100	20	33	<b>FN-510</b>	<b>5.10</b>
150	28	33	<b>FN-511</b>	<b>6.40</b>



### PLASTIC FUNNELS WITH TAPERED STEM

The first three sizes listed are autoclavable. The FN-542 is not.

Diam., Top (mm)	Outside Diam., Stem (mm)	Length, Stem (mm)	NCL #	Price/ea
65	10	21	<b>FN-538</b>	<b>2.80</b>
84	10	23	<b>FN-539</b>	<b>3.30</b>
103	11	34	<b>FN-540</b>	<b>3.40</b>
122	13	44	<b>FN-542</b>	<b>6.45</b>



## PLASTIC FUNNELS WITH NARROW, STANDARD-LENGTH STEM



Top Diameter	Outside Stem Diameter	Stem Length	NCL #	Price/ea
34 mm	6 mm	52 mm	<b>FN-971</b>	<b>2.70</b>
48 mm	7 mm	50 mm	<b>FN-972</b>	<b>2.85</b>
55 mm	7 mm	61 mm	<b>FN-973</b>	<b>2.95</b>
66 mm	7 mm	65 mm	<b>FN-974</b>	<b>3.50</b>
77 mm	7 mm	80 mm	<b>FN-975</b>	<b>4.90</b>
92 mm	8 mm	88 mm	<b>FN-976</b>	<b>5.25</b>
104 mm	9 mm	89 mm	<b>FN-977</b>	<b>7.00</b>

## HIGH SPEED PLASTIC FUNNELS

Unsurpassed for non-vacuum, high-speed filtration, these clear funnels have internal helicoid ribs to increase the filtration speed. Vertical outer ribs prevent airlock. Made from autoclavable polymethylpentene, the funnels are unaffected by corrosive liquids and are heat resistant to 177°C (350°F).



Top Diameter	Outside Stem Diameter	Stem Length	NCL #	Price/ea
51 mm	7.5 mm	150 mm	<b>FB-661</b>	<b>3.85</b>
70 mm	7.5 mm	150 mm	<b>FB-662</b>	<b>4.35</b>
100 mm	10 mm	108 mm	<b>FB-663</b>	<b>6.10</b>
140 mm	16 mm	132 mm	<b>FB-664</b>	<b>21.40</b>

## LARGE FUNNELS

Ideal for filling drums, barrels, tanks, and carboys. Provides ample volume and fast flow. The two smaller sizes have a handle at the rim. Molded natural polypropylene funnels are autoclavable at 121°C (250°F).

Top Diameter	Outside Stem Diameter	Height	NCL #	Price/ea
7.9 in	0.87 in	7.9 in	<b>FB-410</b>	<b>14.00</b>
9.8 in	1.18 in	10.2 in	<b>FB-411</b>	<b>24.00</b>
13.8 in	1.37 in	17.3 in	<b>FB-412</b>	<b>54.65</b>





**GLASS BURET FUNNEL**  
FK-500      8.10/ea



**PLASTIC BURET FUNNEL**  
BN-971      2.95/ea



**DISPOSABLE FILTER FUNNELS WITH RECEIVER BASE**

Sterile, individually wrapped, 150 ml graduated filter funnel with receiver. Includes 47 mm membrane filter. After filtering the sample, upper section twists off for easy removal of membrane filter. For use in coliform determinations where autoclaving is not an option.

Pore Size	Membrane	NCL #	Price/cs
0.2 µm	Plain	<b>FN-120</b>	<b>490.80/72</b>
0.45 µm	Grided	<b>FN-130</b>	<b>490.80/72</b>



**DISPOSABLE FILTER FUNNELS**

Sterile, individually wrapped graduated funnel with 47 mm, 0.45 µm, grided membrane filter. Each case includes 50 funnels and 6 funnel adapters (for use with NCL# R-777, rubber stopper). After filtering the sample, upper section twists off for easy removal of membrane filter.

Funnel Size	NCL #	Price/cs
100 ml	<b>FN-145A</b>	<b>208.00</b>
250 ml	<b>FN-145B</b>	<b>229.20</b>

# FILTER FUNNELS



## MAGNETIC FILTER FUNNEL

Unbreakable polysulfone 300 ml funnel body accommodates 47 mm (4.7 cm) filters. Comes with rubber stopper to fit 1000 or 2000 ml suction flasks, or can be used with a filtration manifold.

**FG-710      263.15/ea**



## KIMBLE GLASS FILTER FUNNEL ASSEMBLY

Supplied with a coarse fritted glass support base, 300 ml graduated funnel, anodized aluminum clamp, and #8 silicone stopper.

**FK-937      154.15/ea**

REPLACEMENT PARTS FOR KIMBLE FILTER ASSEMBLY		
Description	NCL #	Price/ea
glass funnel, 47 mm, 300 ml	<b>FK-937A</b>	<b>63.00</b>
glass support base	<b>FK-937B</b>	<b>89.70</b>
clamp	<b>FK-937C</b>	<b>74.65</b>
silicone stopper with hole	<b>FK-937D</b>	<b>18.75</b>



## MFS GLASS FILTER FUNNEL ASSEMBLY

Supplied with stainless steel screen support base, 300 ml graduated funnel, anodized aluminum clamp, and #8 silicone stopper.

**FS-315      206.00/ea**

## MFS POLYSULFONE FILTER FUNNELS

300 ml graduated polysulfone funnel with recessed filter support to provide sufficient clearance so it will not twist or tear membrane filter. Comes with #8 stopper to fit 1000 ml suction flasks or MFS vacuum manifold. For use with 47 mm filters.

Funnel top locks onto base with  $\frac{1}{4}$  turn.

Probably not as nice as the magnetic version, but a lot cheaper!

**FS-5010      55.70/ea**

