

### GLASS TUBING CUTTER

With brass compression spring which holds jaws open when not in use. Slotted screw retains easily replaceable cutting wheel. For cutting glass tubing up to 30 mm in diameter and cuts from the outside.

**GH-100**                      **15.00/ea**



### BENCH MOUNTED HEAVY DUTY GLASS TUBING CUTTER

Attaches to laboratory countertops up to 1<sup>3</sup>/<sub>4</sub>" thick. Cuts glass tubing up to 1<sup>1</sup>/<sub>2</sub>" in diameter. Aluminum alloy casting.

**GH-250**                      **181.60/ea**

### REPLACEMENT WHEELS FOR GLASS CUTTERS

Packs of 10 wheels.

**GH-101**                      **10.80/pk**

### CARBIDE REPLACEMENT WHEELS FOR GLASS CUTTERS

Packs of 10 wheels.

**GH-102**                      **77.85/pk**



### BROKEN GLASS DISPOSAL BOX

Made of sturdy corrugated cardboard, these receptacles are supplied with a 2 mil thick polyethylene bag to contain broken glass fragments and spilled liquids. A pair of tabs hold the liner in place. When box is full, the safety cap is pulled into place and the entire unit is ready for disposal.

**DB-300**                      8" W x 8" D x 10" H                      **50.75/pk of 6**

**DB-301**                      12" W x 12" D x 27" H                      **67.25/pk of 6**



### GLASS WOOL

Made of soft spun glass. Chemically inert, with no resins or binders. Comes in 1 lb bags.

**GW-200**                      **73.00/ea**