

CIMAREC HOT PLATES



- New LED digital temperature display
- Hot surface warning system
- Low-profile design
- Acid, alkali-resistant seamless ceramic top
- High wattage heating element
- Integral ring stand holder
- Detachable cord
- Microprocessor control offers more consistent slow-speed stirring *
- Stir speed adjustable from 50 to 1200 rpm *
- For 120VAC operation

* applies to stirring hot plates only

STIRRING HOT PLATES

Watts	Size (L x W x H)	NCL #	Price/ea
410	4.25" x 4.25" x 3.6"	HT-015	403.75
1040	7.25" x 7.25" x 3.8"	HT-325	432.25
1440	10.25" x 10.25" x 4"	HT-635	568.10

HOT PLATES ONLY

Watts	Size (L x W x H)	NCL #	Price/ea
410	4.25" x 4.25" x 3.6"	HT-915	226.10
1040	7.25" x 7.25" x 3.8"	HT-225	247.00
1440	10.25" x 10.25" x 4"	HT-535	421.80

STIR PLATES & HOT PLATES WITH ALUMINUM TOP



ALUMINUM-TOP STIRRER, HOT PLATE, AND COMBO

- Rugged 7.25" x 7.25" aluminum surface provides uniform heat distribution and added durability.
- Microprocessor maintains consistent, repeatable temperature settings, ideal for repetitive procedures.
- Easy-to-read digital display indicates temperature settings.
- Bright red warning lights signal when heating surface is above 50°C.
- StirTrac technology offers improved slow-speed stirring, consistent speed control and stronger magnetic coupling.
- Integral ring-stand holder accommodates 0.5" (1.3 cm) diameter support rod.
- All models supplied with detachable cord. 115 volts.

Item	NCL#	Price/ea
Stir Plate	ST-125A	337.25
Stirring hot plate	HT-325A	519.65
Hot plate	HT-225A	319.20



NUOVA II HOT PLATE/STIRRER, (MODEL #18400)

49 sq. in. of white, glass coated stainless steel working space for high chemical resistance and easy cleaning. Continuous stirring speed control from 100 to 1000 rpm. Heat range from 100 to 700°F. Indicator lights for both stirring and heating modes.

HT-184

627.95/ea

CORNING HOT PLATES



- Digital displays indicate when set temperature is reached
- Power button warns if unit is plugged in
- Hot top indicator indicates when top is too hot to touch
- Temperature accuracy is within $\pm 2^{\circ}\text{C}$ with use of optional temperature controlling probe

Heating Surface	NCL #	Price/ea
5" x 7"	PC-400D	236.90
10" x 10"	PC-600D	405.40



CORNING STIRRING HOT PLATES

- Dual digital displays show temperature and stir speed
- Power button warns if unit is plugged in
- Hot top indicator indicates when top is too hot to touch
- Temperature accuracy is within $\pm 2^{\circ}\text{C}$ with use of optional temperature controlling probe

Heating Surface	NCL #	Price/ea
5" x 7"	PC-420D	415.90
10" x 10"	PC-620D	547.50

Optional temperature control probe for Corning digital hot plates and hot plate stirrers.
PC-679 **123.40/ea**