


## REAGENT BOTTLE LABELS




**DATE RECEIVED** \_\_\_\_\_  
**DATE OPENED** \_\_\_\_\_  
 Reorder 1.800.648.7836

For use on reagents which you purchase. Label has space for user to write date received and date opened. Label size is  $\frac{3}{4}$ " high by  $1\frac{1}{2}$ " long. Die cut and self-adhesive. 500 labels per roll.

**RL-100 9.00/rl**



**LABORATORY REAGENT**  
 PREP. DATE \_\_\_\_\_ CHECKED \_\_\_\_\_  
 EXPIRATION DATE \_\_\_\_\_  
 REAGENT \_\_\_\_\_  
 CONC. \_\_\_\_\_  
 PROCEDURE \_\_\_\_\_  
 MADE BY \_\_\_\_\_

4" wide x 2" high designed for use of reagents made in your lab. Has space for "Prep. date," "checked," "exp. date," "reagent," "conc.," "procedure," and "made by". Self-adhesive, 500 per roll.


**RL-150 35.70/rl**



## BIOHAZARD LABEL (OSHA)

Permanent, pressure-sensitive, self-adhesive Biohazard label measures  $2\frac{1}{2}$ " wide x 3" high. Brilliant fluorescent color as mandated by OSHA. 500 per roll.

**HL-100 27.00/rl**



## SINGLE-ROLL TAPE DISPENSER

Accommodates tape rolls from  $\frac{1}{2}$ " wide to 1" wide with 1" cores.

**LB-111 3.85/ea**



## MICROSCOPE SLIDE LABELS

Labels fit 1 x 3" (25 x 75 mm) standard microscope slide. Labels are  $\frac{15}{16}$ " x  $\frac{15}{16}$ " with slightly rounded corners. White, die-cut, self adhesive, 100 per roll.


**LS-150 19.75/rl**

## HFR WARNING LABELS

<b>H</b> HEALTH	<input type="checkbox"/>	OWNER
<b>F</b> FLAMMABILITY	<input type="checkbox"/>	LOT NO.
<b>R</b> REACTIVITY	<input type="checkbox"/>	REC'D.
PERSONAL PROTECTION	<input type="checkbox"/>	EXP.

2" x 2" as pictured above. Self-adhesive, die-cut, 500 per roll.

**HL-160 27.00/rl**



Chemical Name \_\_\_\_\_  
 Common Name \_\_\_\_\_  
 Manufacturer \_\_\_\_\_  
 MSDS# \_\_\_\_\_ Date \_\_\_\_\_

1" x 3" as pictured above. Self-adhesive, die-cut, 500 per roll.

**HL-170 22.80/rl**


## LABELING TAPE DISPENSER, WRITE ON

This unique tape dispenser features a convenient writing platform which allows both right-handed and left-handed users to write on the tape before applying it to the item. The epoxy coated platform has finger holes to facilitate peeling and a built-in cutting edge. Has a weighted base and non-skid bottom. Includes a plastic core for 3" diameter rolls of tape up to 1" wide. Tape not included.

**LB-100 71.70/ea**



Labeling Tape			
Autoclavable, self-adhesive. Rolls are 40 yards (1440") long.			
Width	Color	NCL #	Price/roll
1/2"	white	<b>LV-102A</b>	<b>9.40</b>
3/4"	white	<b>LV-103A</b>	<b>12.50</b>
1"	white	<b>LV-104A</b>	<b>15.60</b>
1/2"	blue	<b>LV-102B</b>	<b>9.40</b>
3/4"	blue	<b>LV-103B</b>	<b>12.50</b>
1"	blue	<b>LV-104B</b>	<b>15.60</b>
1/2"	red	<b>LV-102R</b>	<b>9.40</b>
3/4"	red	<b>LV-103R</b>	<b>12.50</b>
1"	red	<b>LV-104R</b>	<b>15.60</b>
1/2"	yellow	<b>LV-102Y</b>	<b>9.40</b>
3/4"	yellow	<b>LV-103Y</b>	<b>12.50</b>
1"	yellow	<b>LV-104Y</b>	<b>15.60</b>



Labels, Self-Adhesive			
Individually pre-cut self-adhesive labels with red borders.			
Size	Quantity/roll	NCL #	Price/roll
1" x 3/4"	500	<b>LF-100</b>	<b>9.20</b>
1" x 1 5/8"	500	<b>LF-101</b>	<b>12.15</b>
1 1/2" x 2 1/2"	500	<b>LF-102</b>	<b>15.75</b>

## WRITE-ON ALUMINUM TAGS

Permanently engrave your messages on metal tags with an ordinary ballpoint pen. Extremely durable and suited for use in hot, cold, and moist environments. Handy size of 3 1/2" x 1 3/8" makes it ideal for tagging samples, cylinders, valves, tubing, etc.

Ideal for tagging your thermometers with number and correction factor after comparing them to a certified thermometer.

**TC-307 6.20/pk of 6 40.00/cs of 48**

