


# HERATHERM OVENS



Thermo Scientific newly designed laboratory ovens. Feature advanced digital controllers for optimum temperature precision. Unit with mechanical convection has better uniformity. Supplied with 2 shelves, holds a maximum of 13 shelves. Supplied with power cord.

<b>Model #</b>	OMH60	OGH60
<b>Air Circulation</b>	Mechanical	Gravity
<b>Temperature Range</b>	50°C to 330°C	50°C to 330°C
<b>Inside Dimensions</b>	13.9" W x 13.5" D x 20" H	12.9" W x 15.3" D x 18.9" H
<b>Outside Dimensions</b>	20.9" W x 22.2" D x 28.3" H	20.9" W x 22.2" D x 28.3" H
<b>Sensitivity</b>	±0.2°C	±0.3°C
<b>Uniformity</b>	±1.5°C	±2.5°C
<b>Controller</b>	microprocessor with LED display of set-point & temp	microprocessor with LED display of set-point & temp
<b>Electrical</b>	120V, 1440 watts, 12 amps	120V, 1730 watts, 14.4 amps
<b>Shipping Weight</b>	97 lbs.	97 lbs.
<b>NCL #</b>	<b>OV-124M</b>	<b>OV-115G</b>
<b>Price</b>	<b>2425.50/ea</b>	<b>1695.00/ea</b>