

PETRI DISHES (DISPOSABLE)

Note: All petri dishes listed on this page are approximately 50 mm in diameter, and are designed to hold 47 mm diameter filters and pads. Another mystery solved.

Most people consider petri dishes to be a "generic" commodity. That is, all petri dishes are pretty much equivalent, and choices between various brands are a matter of personal preference.

There are many packaging choices with petri dishes, grid filters, and pads (You can get petri dishes with or without pads, grid filters with or without pads, or all three separately). Our recommendation is to get the petri dishes preloaded with pads, and to get the grid filters individually wrapped and sterile. This usually turns out to be the cheapest alternative, plus, chances of accidental contamination are reduced because you don't have to handle the pad. Again, this is a matter of personal preference.

PETRI DISHES WITH PADS, 50 MM Sterile dishes pre-loaded with 47 mm pads.		
Manufacturer	NCL#	Price
Millipore	PM-150	65.55/pk of 150
MFS	PS-100	33.95/pk of 100
Millipore	PM-150A	197.60/cs of 600
MFS	PS-100A	152.65/cs of 500
PETRI DISHES WITHOUT PADS, 50 MM Sterile dishes - NO pads.		
Manufacturer	NCL#	Price
Millipore	PM-151	51.30/pk of 150
MFS	PS-101	29.85/pk of 100
Millipore	PM-151A	151.05/cs of 600
MFS	PS-101A	134.40/cs of 500



ABSORBENT PAD KIT Pads only for 47 to 50 mm petri dishes.		
Description	NCL #	Price
Case of 1000 with dispenser	FG-700	170.65/cs
Pack of 100 (no dispenser)	FG-700B	21.00/pk
Dispenser only	FG-700C	3.00/ea

DISPOSABLE PETRI DISHES

100x15 mm sterile polystyrene stackable petri dishes with triple vented lid. Recommended for microbiology use.

P-200 89.00/cs of 500

