

ACCESSORIES FOR **THERMO** SCIENTIFIC ORION STAR SERIES METERS

THERMO SCIENTIFIC ORION TRIODES (pH/ATC) PROBES FOR STAR SERIES METERS ONLY

Epoxy body for routine pH measurement. Probes have 3 ft cables, except for the waterproof probe, which has a 10 ft cable.

Probe Style	NCL #	Price/ea
gel-filled standard	O-9107BNMD	212.40
refillable standard	O-9157BNMD	297.90
waterproof, gel filled	O-9107WMMD	223.20

THERMO SCIENTIFIC ORION ATC PROBES FOR STAR SERIES METERS ONLY

Not compatible with other series meters.

Body	NCL #	Price/ea
epoxy	O-927005MD	224.10
glass	O-927006MD	239.40
stainless steel	O-927007MD	270.00

pH ACCESSORIES FOR OLDER ORION METERS

The following accessories are compatible only with **Thermo** Scientific Orion "A" and "A+" series meters.

THERMO SCIENTIFIC ORION ATC PROBES

The following ATC probes will work with Orion Models 210A, 230A, 290A, 420A, 720A, and 920A or **Thermo** Scientific Orion Models 210A+, 230A+, 290A+, 420A+, 720A+, and 920A+.

Body	NCL #	Price/ea
epoxy	O-917005	204.30
glass	O-917006	204.30
stainless steel	O-917007	239.40

THERMO SCIENTIFIC ORION TRIODES (pH/ATC) PROBES WITH BNC CONNECTORS

Probe Style	Body	Recommended Use	NCL #	Price/ea
gel-filled standard	epoxy	routine pH	O-9107BN	198.00
refillable standard	epoxy	routine pH	O-9157BN	292.50