

PROBES FOR THE YSI PRO20

Description	NCL #	Price/ea
Self-stirring BOD probe with 2 m (6 ft) cable included	YP-5780	660.25
Galvanic Sensor for any cable	YP-5202	156.75
Polarographic Sensor for any cable	YP-5203	156.75

PROBES AND MEMBRANES FOR THE YSI-200

Description	NCL #	Price/ea
DO/Temp probe with 12 ft. cable	DY-2004	251.75
DO/Temp probe with 30 ft. cable	DY-2010	280.25
Six screw-on membrane caps & 1 bottle of fill solution	YSI-5908	46.80

Membrane kits				
For	Contains	NCL #	YSI #	Price/kit
5905, 5010, 550, 5239, 556, 85	6 threaded membrane caps, 1 btl fill sol'n (black caps)	YSI-5906	5906	46.80
550A, DO-200, 556 (w/older software)	6 threaded membrane caps, 1 btl fill sol'n (yellow caps)	YSI-5908	5908	46.80
556 (w/1.03 software)	6 threaded membrane caps, 1 btl fill sol'n (blue caps)	YSI-5909	5909	48.60
55, 5739, 5750	30 reg membranes, 1 btl filling sol'n	DY-850	5775	29.50
55, 5739, 5750	30 reg membranes, 1 btl filling sol'n	DY-851	5776	37.00

Reconditioning kits				
For	Contains	NCL #	YSI #	Price/kit
55, 5739, 5750	10 sanding disks, sanding tool	DY-858	5680	30.00
550. 550A	10 sanding disks, sanding tool	YSI-5238	5238	29.00



**SELF-STIRRING BOD PROBE,
YSI MODEL 5905**

Features threaded membrane caps and solid stirring rod, eliminating the need to play around with loose membranes, O-rings, and beater boots. Compatible with YSI 51B, 54ARC, 50, 50B, 52, 57, 58, and 59 meters. Requires optional YSI-5011 adapter for compatibility with Models 5000 and 5100. Includes cable, six membrane caps (YSI-5906), and a bottle of filling solution.

YSI-5905 595.00/ea

**SELF-STIRRING BOD PROBE,
YSI MODEL 5010**

Same as the above 5905 probe, except that cable fits YSI models 5000 and 5100 DO meters. Not for use with any other DO meters.

YSI-5010 562.50/ea

**CABLE ADAPTER FOR ADAPTING YSI
5905 PROBE TO 5000 & 5100 DO METERS**

YSI-5011 90.00/ea



**NON-STIRRING BOD PROBE,
YSI 5750**

Fits standard 300 ml BOD bottles. Requires sample agitation with magnetic stirrer for correct operation. Comes with a 5-foot cable and one 5775 membrane kit. Requires optional YSI-5011 adapter for compatibility with Models 5000 and 5100.

DY-835 356.25/ea

**YSI FIELD PROBE WITH SCREW-
ON MEMBRANE CAPS**



New YSI field probe features screw-on membrane caps instead of loose membranes. All probes have cable attached. Choice of 10, 25, or 50-foot cables. Compatible with YSI Model 50, 50B, 51B, 52, 54, 57, and 58 meters. Include one pack membrane caps (YSI-5906).

Cable Length	NCL #	Price/ea
10 feet	YSI-5239A	445.50
25 feet	YSI-5239B	472.50
50 feet	YSI-5239C	517.50



FIELD PROBE, YSI 5739

The standard probe for field DO measurements. Easily removable guard protects the sensor tip. Choose cable length you desire and purchase separately.

DY-832 308.75/ea

DO PROBE CABLES

For use with YSI 5739 Field Probe. Cables longer than 10 ft include reels.

Cable Length	NCL #	Price/ea
10 feet	DY-840	237.50
25 feet	DY-841	261.25
50 feet	DY-842	308.75



ORION® 3-STAR BENCHTOP DO

BENCHTOP DO KIT

Includes meter, 081010MD probe with BOD funnel, 080113 probe service kit, electrode stand, 080017 DO cal sleeve.

DO-1113A **1546.20/ea**

BOD KIT

Includes meter, 086030MD autostir DO probe, 080514 probe electrolyte solution.

DO-1114 **1678.50/ea**

NON-STIRRING REPLACEMENT PROBE

With BOD funnel. For Star series DO meters.

O-081010MD **436.50/ea**

SELF-STIRRING REPLACEMENT PROBE

Funnel is built into probe. For Star series DO meters.

O-086030MD **1283.40/ea**

ORION® 3-STAR PORTABLE DO

PORTABLE DO KIT

Includes meter, 083010MD epoxy field probe with 10 ft cable, 080513 DO probe maintenance kit, 080017 DO cal sleeve and hard field case.

DO-1213A **1187.10/ea**

ORION® 4-STAR PORTABLE pH/DO

PORTABLE pH/DO METER KIT II

Same as above, except different probes (with 4 ft cables): 083005MD DO probe and 9017BNMD pH/ATC triode.

DP-1216A **1339.20/ea**

3-STAR RDO METER (RDO ONLY)

Portable Meter only

RO-1213 **535.50/ea**

Portable Meter, RDO probe with 3 m cable

RO-1213A **1106.10/ea**

Portable Meter, RDO probe with 6 m cable

RO-1213B **1229.40/ea**

4-STAR RDO METER (RDO AND pH)

Portable Meter only

PR-1213 **868.50/ea**

Portable Meter, RDO probe with 3 m cable

PR-1213A **1229.40/ea**

Portable Meter, RDO probe with 6 m cable

PR-1213B **1376.10/ea**

5-STAR RDO METER (RDO, pH, AND CONDUCTIVITY)

Meter only

MO-1213 **963.90/ea**

Meter, RDO probe with 3 m cable

MO-1213A **2162.70/ea**

Meter, RDO probe with 6 m cable

MO-1213B **2410.20/ea**

PROBE ONLY

with 3 m cable **O-087010** **646.20/ea**

with 6 m cable **O-087020** **695.70/ea**

with 10 m cable **O-087030** **721.80/ea**

with 15 m cable **O-087050** **746.10/ea**

with 30 m cable **O-087100** **794.70/ea**

REPLACEMENT OPTICAL CAP

O-087001 **80.10/ea**

REPLACEMENT PROBE GUARD

O-087002 **296.10/ea**

CALIBRATION SLEEVE

O-087003 **114.30/ea**



THERMO SCIENTIFIC ORION DO PROBE

Determine the DO using any good pH meter. Simple to use, accurate, and no more messing around with loose membranes or filling solutions. Built-in magnetic stir bar. Requires magnetic stir plate for correct operation. Will work with Thermo Scientific Orion, Corning, or any other good pH meter, but not compatible with YSI meters. For laboratory use only, not for field use. Includes two bonded membrane caps which usually last 3 to 6 months each.

Replacement Parts for Thermo Scientific Orion DO Probe			
Description	Orion #	NCL #	Price/ea
Bonded Membrane Caps	970801	DO-892	63.00
Batteries	970803	DO-893	8.00
Funnel	970802	DO-894	134.90

Connector Type	BNC
NCL #	DO-891
Orion #	970899
Price	1085.40/ea

REPLACEMENT PROBES FOR ORION STAR SERIES DO METERS

THERMO SCIENTIFIC ORION DO PROBE
O-081010MD 436.50/ea



THERMO SCIENTIFIC ORION
FIELD DO PROBE
O-083010MD 665.00/ea



THERMO SCIENTIFIC
ORION SELF-STIRRING
BOD PROBE
O-086030MD 1283.40/ea



Description	For DO Probes...	Orion #	NCL #	Price/ea
DO Maintenance Kit A	O-083010MD, O-086030MD	080513	O-080513	144.90
1 Membrane Cap		080515	O-080515	40.50
DO Maintenance Kit B	O-081010MD	080113	O-080113	91.80
3 Membrane Caps		081003	O-081003	91.80



YSI OPTICAL DO METER ACCESSORIES



OPTICAL FIELD PROBE CAP

Replacement cap for YSI Optical Field Probe.
1 year warranty.*

*warranty applies to manufacturing defects only.

YSI-320 91.00/ea




OPTICAL BOD PROBE CAP

Replacement cap for YSI Optical BOD Probe.
Not for use with YSI Optical Field Probe.

1 year warranty.*


YSI-482 95.00/ea



PROBOD OPTICAL SELF-STIRRING PROBE

Smart Sensor retains calibration in the probe. Does not require warm up or electrode maintenance. Replaceable stir rod. Sensor cap included with initial probe purchase. Two year warranty on probe and cable.

YSI-400 850.25/ea



YSI BOD ANALYST PRO SOFTWARE

Connect meter directly to the PC and BOD Analyst Pro software to take initial and final DO values and let the software eliminate human errors by calculating the BOD/CBOD results. Data can also be entered manually without connection to a computer. For use with YSI instruments, 52, 5000, 5100, Professional Plus, or ProODO.

YSI-120 495.00/ea



METER DOCK WITH BOD BOTTLE HOLDER

For use with all YSI Pro Series Meters.

YSI-444 44.00/ea



USB POWER SUPPLY

For use with all YSI Pro Series Meters with USB saddle.

YSI-518 40.00/ea



TRIPOD METER STAND

For use with all YSI Pro Series meters.

YSI-507 26.00/ea