

THERMO SCIENTIFIC ORION FLUORIDE IONPLUS COMBINATION PROBE

For rapid determinations of Fluoride in water. Polyphosphates do not interfere. No separate reference probe required. Includes probe and one 60 ml bottle of filling solution. BNC connector.

O-9609BN 668.70/ea

FILLING SOLUTION FOR O-9609BN ONLY

O-900061 18.80/ea 84.60/pk of 5x60 ml



THERMO SCIENTIFIC ORION FLUORIDE PROBE. BNC connector. Requires 0-900100

Description	NCL #	Price/ea
Fluoride probe	O-9409BN	564.30
Single junction reference probe	O-900100	178.20

Filling Solution for O-900100 when used with O-9409BN Fluoride Probe **ONLY**
O-900001 17.00/ea 76.50/pk of 5x60 ml

ISA Buffers for use with Fluoride Probes

Please see the chemical section in this catalog for Fluoride buffers in other sizes.

Concentration	Manufactured by	NCL #	Price
TISAB III Concentrate	Thermo Scientific Orion	O-940911	121.50/475 ml HM
TISAB III Concentrate	NCL	F-29	33.00/500 ml HM
TISAB II	Thermo Scientific Orion	O-940909	112.50/3.8 liter
TISAB II	NCL	F-29A	55.00/3.8 liter

Note: TISAB II is mixed 1:1 with samples and standards. TISAB III is mixed 1:10 with samples and standards. Therefore, TISAB III would last you 10 times as long as TISAB II.