


ORP PROBE

THERMO SCIENTIFIC ORION ORP/REDOX PROBE

For measuring Oxidation-Reduction Potential. No reference probe required. Works with any pH meter which has a millivolt (mv) mode.



Description	NCL #	Price
Sure-Flow with BNC connector	O-9678BN	393.30/ea
Gel filled triode for A"+"&" meters	O-9179BN	279.00/ea
Refillable triode for A"+"&" meters	O-9180BN	365.00/ea
Gel filled triode for Star Meters	O-9179BNMD	305.00/ea
Refillable triode for Star Meters	O-9180BNMD	389.00/ea



PORTABLE WATERPROOF ORP PACKAGE WITH SUBMERSIBLE ORP PROBE

Portable waterproof meter developed for hard on-site work. Features a 200-point datalogger, large LCD display, and a RS232 interface. Battery operated using four AA disposable batteries. Datalogging intervals are between 5 seconds and one hour. Logged data can be viewed from meter display or can be downloaded to a PC using optional software and cable. Package contains meter, case, and a submersible ORP probe with 20-foot cable. 3-Year warranty on meter.

Note: The filling solution for the probe is saturated KCl (approx. 3.8 M).

Description	NCL #	Price/ea
Meter, Submersible ORP Probe with 20 ft cable, and case	NCL-100	985.00
Replacement ORP Probe with 20 ft cable	NCL-101	228.80
Data collection software and cable	WT-852	175.00