

SOLID RUBBER STOPPERS (NO HOLES)

Solid black rubber. Order individually or by the pack.
Each pack contains approximately one pound.



Size (#)	Top (mm)	Bottom (mm)	Height (mm)	Quantity/pk	NCL#	Price/ea	Price/pk
000	13	8	21	180	R-000	.14	22.00
00	15	10	25	117	R-00	.24	22.00
0	17	13	25	82	R-0	.32	22.00
1	19	14	25	57	R-1	.44	22.00
2	20	16	25	52	R-2	.48	22.00
3	24	18	25	39	R-3	.64	22.00
4	26	20	25	33	R-4	.76	22.00
5	27	23	25	28	R-5	.89	22.00
5.5	28	24	25	26	R-5.5	.96	22.00
6	32	26	25	20	R-6	1.25	22.00
6.5	34	27	25	18	R-6.5	1.40	22.00
7	37	30	25	15	R-7	1.66	22.00
7.5	39	31	25	14	R-7.5	1.93	22.00
8	41	33	25	12	R-8	2.08	22.00
8.5	43	36	25	11	R-8.5	2.38	22.00
9	45	37	25	10	R-9	2.50	22.00
9.5	46	38	25	9	R-9.5	2.88	22.00
10	50	42	25	7	R-10	3.56	22.00
10.5	53	45	25	7	R-10.5	3.73	22.00
11	56	48	26	6	R-11	4.15	22.00
11.5	63	50	25	5	R-11.5	5.14	22.00
12	64	54	25	4	R-12	6.25	22.00
15	103	83	39	1	R-15	18.15	N/A

ONE HOLE RUBBER STOPPERS

Holes are approximately 3 mm in diameter. If you need larger holes, you will need to drill them out. Order individually or by the pack. Each pack contains approximately one pound.



Size (#)	Top (mm)	Bottom (mm)	Height (mm)	Quantity/pk	NCL#	Price/ea	Price/pk
000	13	8	21	180	R-000A	.14	22.00
00	15	10	25	117	R-00A	.32	22.00
0	17	13	25	82	R-0A	.32	22.00
1	19	14	25	57	R-1A	.50	22.00
2	20	16	25	52	R-2A	.53	22.00
3	24	18	25	39	R-3A	.68	22.00
4	26	20	25	33	R-4A	.81	22.00
5	27	23	25	28	R-5A	.98	22.00
5.5	28	24	25	26	R-5.5A	.96	22.00
6	32	26	25	20	R-6A	1.35	22.00
6.5	34	27	25	18	R-6.5A	1.50	22.00
7	37	30	25	15	R-7A	1.93	22.00
7.5	39	31	25	14	R-7.5A	1.93	22.00
8	41	33	25	12	R-8A	2.38	22.00
8.5	43	36	25	11	R-8.5A	2.38	22.00
9	45	37	25	10	R-9A	2.86	22.00
9.5	46	38	25	9	R-9.5A	2.92	22.00
10	50	42	25	7	R-10A	3.56	22.00
10.5	53	45	25	7	R-10.5A	4.18	22.00
11	56	48	26	6	R-11A	4.90	22.00
11.5	63	50	25	5	R-11.5A	5.16	22.00
12	64	54	25	4	R-12A	6.25	22.00

2-HOLE RUBBER STOPPERS

Holes are approximately 3 mm in diameter. If you need larger holes, you will need to drill them out. Order individually or by the pack. Each pack contains approximately one pound.



Size (#)	Top (mm)	Bottom (mm)	Height (mm)	Quantity/pk	NCL#	Price/ea	Price/pk
000	13	8	21	180	R-000B	.14	22.00
00	15	10	25	117	R-00B	.32	22.00
0	17	13	25	82	R-0B	.32	22.00
1	19	14	25	57	R-1B	.50	22.00
2	20	16	25	52	R-2B	.53	22.00
3	24	18	25	39	R-3B	.68	22.00
4	26	20	25	33	R-4B	.81	22.00
5	27	23	25	28	R-5B	.98	22.00
5.5	28	24	25	26	R-5.5B	.96	22.00
6	32	26	25	20	R-6B	1.35	22.00
6.5	34	27	25	18	R-6.5B	1.50	22.00
7	37	30	25	15	R-7B	1.93	22.00
7.5	39	31	25	14	R-7.5B	1.93	22.00
8	41	33	25	12	R-8B	2.38	22.00
8.5	43	36	25	11	R-8.5B	2.38	22.00
9	45	37	25	10	R-9B	2.86	22.00
9.5	46	38	25	9	R-9.5B	2.92	22.00
10	50	42	25	7	R-10B	3.56	22.00
10.5	53	45	25	7	R-10.5B	4.18	22.00
11	56	48	26	6	R-11B	4.90	22.00
11.5	63	50	25	5	R-11.5B	5.16	22.00
12	64	54	25	4	R-12B	6.25	22.00

STOPPER SIZER

Two scale stopper sizer determines the correct size of the rubber (sizes 00 to 12) or cork (sizes 0 to 26) stopper needed for any flask, bottle, etc. Just insert the sizer into the vessel and the correct stopper size is read at the top of the neck. For openings up to 66 mm (2³/₈").

RS-202 9.00/ea



LARGE HOLE RUBBER STOPPER

Size 8 rubber stopper has single hole approx. 1/2" (12 mm) in diameter which allows it to fit many Buchner and filter funnels.

R-777 9.26/ea 20.85/pk of 3

CORK BORER SET

Brass plated cork borer set, consisting of 15 sizes from 4.5 mm to 25 mm. Punch included. Borers nest for easy storage.

CF-845 49.00/set

SHARPENER FOR ABOVE

CF-888 39.00/ea

