

Lightweight, foam bottle shippers for safe and economical transport of single or multi-bottle shipments. Molded two-piece units hold glass bottles securely and eliminate breakage due to shock and vibration. All have been laboratory tested and accepted by UPS, Parcel Post, and common carriers.

Fast assembly with interlocking edges and molded-in tape grooves. Pressure-sensitive labels (not included) adhere directly to foam shippers. No outer cartons are necessary. Bottles not included.



CP-338

STYROFOAM SHIPPERS



CP-393

For Bottle Size	Inside Dimensions	Outside Dimensions	NCL#	Price/ea	Price/cs
1 Gallon (Cider Jug)	6 ¹ / ₄ " dia. x 12 ³ / ₈ " H	10" x 10" x 15 ¹ / ₈ "	CP-338	9.20	58.90/8
1 Gallon (Pharmaceutical)	6 ¹ / ₄ " dia. x 13 ⁵ / ₈ " H	9 ¹ / ₄ " x 9 ¹ / ₄ " x 16 ¹ / ₂ "	CP-380	8.60	68.85/9
500 ml	3" dia. x 7 ⁵ / ₈ " H	4 ³ / ₄ " x 4 ³ / ₄ " x 10"	CP-358	2.50	125.10/60
2.5 Liter/80 oz jugs & Boston Rounds	5 ³ / ₈ " dia. x 12 ³ / ₄ " H	6 ⁷ / ₈ " x 6 ⁷ / ₈ " x 14 ¹ / ₄ "	CP-359	4.00	56.95/18
1 Liter/1 Quart Boston Rounds, some wide mouth Quarts	3 ⁷ / ₈ " dia. x 9" H	5 ³ / ₄ " x 5 ¹ / ₂ " x 10 ¹ / ₂ "	CP-393	4.00	145.15/45
1 Pint Boston Rounds	3 ¹ / ₄ " dia. x 7 ¹ / ₂ " H	4 ¹ / ₂ " x 4 ³ / ₄ " x 9 ¹ / ₄ "	CP-361	2.60	140.00/64
Half Pint Boston Round	2 ³ / ₈ " dia. x 5 ⁵ / ₈ " H	4" x 4" x 7"	CP-344	1.65	153.00/120

STYROFOAM SLIDE MAILERS

Slide mailer carries up to ten standard 3" x 1" slides placed back to back in each slot for protection from breakage in transit. Fits CP-370 paperboard mailing sleeve.


CP-369 210.05/cs of 300



MAILING SLEEVES FOR SLIDE MAILERS

Increase the cost effectiveness of these foam mailers by using sturdy mailing sleeves. Strong 18-pt. white paperboard sleeves with reverse-tuck closures enable CP-369 slide mailers to withstand multiple trips.

CP-370 115.15/cs of 300



COLIFORM WATER TEST MAILER

Avoid shipping damage to analytical samples in transit with the newest foam protection. Constructed of 1/2" thick expanded polystyrene, the engineered rib structure is ideal for holding BC-170 sample container firmly in place, insuring sample integrity. A combination of the molded two piece interlocking mailer and the leash on Corning's coliform water test sampling container forces the sample to stay capped during transit. Outside dimensions: 4⁷/₈" x 3³/₄" x 3³/₈". Inside dimensions: 3⁵/₈" x 1¹/₂" x 2¹/₄".

BC-605 99.30/cs of 96

MAILING SLEEVES FOR COLIFORM MAILER

Municipalities, laboratories analyzing water and independent testing agencies will find the mailer essential. It is lightweight, reusable and recyclable. For extra protection, cardboard mailing sleeves are available for the mailer. Sleeves may be printed with laboratory name for a slight additional charge.

BC-606 123.30/cs of 192