



### DRY ACID SPILL KIT

Contains:

- (1) 10 lb Kolor-safe® dry acid neutralizer
- (3) 3" x 4' Hazmat socks
- (10) 12" x 12" Hazmat pads
- (1) scoop with detachable scraper
- (1) pair splash goggles
- (1) pair nitrile gloves
- (2) chemical classifier strips
- (1) chemical classifier chart
- (1) 5 gallon plastic bucket with snap-on lid
- (1) disposal bag & twist tie

**SK-275**

**144.90/ea**



### LIQUID ACID SPILL KIT

Contains:

- (1) gallon Kolor-safe® liquid acid neutralizer
- (3) 3" x 4' Hazmat socks
- (10) 12" x 12" Hazmat pads
- (1) pair splash goggles
- (1) pair nitrile gloves
- (2) chemical classifier strips
- (1) chemical classifier chart
- (1) 5 gallon plastic bucket with snap-on lid
- (1) disposal bag & twist tie

**SK-431**

**99.20/ea**



### DRY BASE SPILL KIT

Contains:

- (1) 10 lb Kolor-safe® dry base neutralizer
- (3) 3" x 4' Hazmat socks
- (10) 12" x 12" Hazmat pads
- (1) scoop with detachable scraper
- (1) pair splash goggles
- (1) pair nitrile gloves
- (2) chemical classifier strips
- (1) chemical classifier chart
- (1) 5 gallon plastic bucket with snap-on lid
- (1) disposal bag & twist tie

**SK-276**

**144.90/ea**



### LIQUID BASE SPILL KIT

Contains:

- (1) gallon Kolor-safe® liquid base neutralizer
- (3) 3" x 4' Hazmat socks
- (10) 12" x 12" Hazmat pads
- (1) pair splash goggles
- (1) pair nitrile gloves
- (2) chemical classifier strips
- (1) chemical classifier chart
- (1) 5 gallon plastic bucket with snap-on lid
- (1) disposal bag & twist tie

**SK-453**

**99.20/ea**

Kits are packed in 5 gallon plastic buckets for easy portability and disposal of spill cleanup materials. Everything needed to quickly stop small spills before they damage property or escape down drains. Features unique Kolor-safe® neutralizers that change color to indicate that neutralization is complete

### PORTABLE SPILL KIT

Designed for portability and quick response to containing all types of non-aggressive liquid spills and leaks. Absorbs up to 4 gallons.

Contains:

- (4) 3" x 4' socks (non-haz)
- (3) 16" x 18" premium pads (non-haz)
- (1) pair nitrile gloves
- (1) pair splash goggles
- (1) 5 gallon plastic bucket

**SK-455**

**62.60/ea**



### HAZMAT SPILL STATION

Designed for larger spills and leaks. Stored in UN-Certified Open-head Poly (with lid).

Contains:

- (12) 3" x 4' Hazmat socks
- (2) 3" x 10' Hazmat socks
- (100) 16" x 18" Hazmat pads
- (12) 8.5" x 17" Hazmat pads
- (1) 10 lb dry acid neutralizer
- (1) 10 lb dry base neutralizer
- (1) 10 lb dry hydrofluoric acid neutralizer
- (1) 1" GapSeal® plugging compound
- (1) pair splash goggles
- (1) pair nitrile gloves
- (1) chemical classifier chart
- (2) chemical classifier strips
- (1) disposal bag & twist tie
- (1) 55 gallon Poly drum with Lever-Lock metal ring top.

**SK-201**

**531.55/ea**



### MERCURY SPILL KIT

Designed to handle all stages of cleaning up mercury spills. The MERCORSORB Spill Kit is designed for one-time use only. MERCORSORB Kits do not contain any chemicals designed to alter the mercury's toxicity or danger. Recommended for use by trained personnel only.

Contains:

- (1) 270 g amalgamation powder
- (1) vapor suppressor shaker bottle
- (1) waste collection bottle
- (1) chemical sponge
- (1) disposal bag & twist tie
- (1) bottle indicator powder (250 g)
- (1) aspirator bottle
- (1) mixing tub with spatula
- (1) pair nitrile gloves
- (1) wisk on pan

**SK-252**

**105.55/ea**



## SPILL KIT REPLACEMENT PARTS

Description	For Kit(s)	Quantity	NCL#	Price
3" x 4' Hazmat socks	SK-201, 275, 276, 431, 453	cs (15/cs)	<b>SR-230</b>	<b>64.35</b>
12" x 12" Hazmat pads	same as above	cs (100/cs)	<b>SR-271</b>	<b>40.50</b>
Dry acid neutralizer, 10#	SK-201, SK-275	each	<b>SR-440</b>	<b>29.70</b>
Dry base neutralizer, 10#	SK-201, SK-276	each	<b>SR-450</b>	<b>46.75</b>
Liquid acid neutralizer	SK-431	cs (4 x 1 gal)	<b>SR-410</b>	<b>88.45</b>
Liquid base neutralizer	SK-453	cs (4 x 1 gal)	<b>SR-420</b>	<b>87.75</b>
Mercury amalgamation powder, 270 grams	SK-252	each	<b>SR-520</b>	<b>15.00</b>
Mercury indicator powder, 250 grams	SK-252	each	<b>SR-523</b>	<b>23.50</b>
Mercury vapor suppressor	SK-252	cs (5 x 1 lb)	<b>SR-460</b>	<b>111.35</b>
3" x 4' non-haz socks	SK-455	cs (15/cs)	<b>SR-130</b>	<b>30.30</b>
16" x 18" non-haz pads	SK-455	cs (100/cs)	<b>SR-175</b>	<b>70.35</b>
16" x 18" Hazmat pads	SK-201	cs (100/cs)	<b>SR-170</b>	<b>112.00</b>
8.5" x 17" Hazmat pillows	SK-201	cs (10/cs)	<b>SR-166</b>	<b>49.50</b>
3" x 10' Hazmat socks	SK-201	cs (6/cs)	<b>SR-137</b>	<b>68.90</b>
Chemical classifier strips	all	cs (10/cs)	<b>SR-570</b>	<b>73.50</b>
Wastewater classifier strips	all	cs (10/cs)	<b>SR-580</b>	<b>73.50</b>

### SPILL TRAY

Eliminate spills with this chemical-resistant, maintenance-free, polyethylene containment tray. Use one at each workstation or fumehood to safely handle chemicals and liquids while pipetting, measuring, or mixing. Holds 1<sup>1</sup>/<sub>4</sub> liters of spilled liquid. Easy to clean—simply lift off the white grid and empty. Low profile design is perfect for drying glassware and delicate items. Plastic grid cushions and protects against breakage and scratches. Grid cuts drying time in half by permitting air circulation. Entire unit is dishwasher safe. Size is 16<sup>1</sup>/<sub>2</sub> x 13 x <sup>5</sup>/<sub>8</sub> inches. Weight is 1 pound.

**CS-345**                      **17.65/ea**

